User manual

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..... PAGE 2

..... PAGE 30

..... PAGE 58

TAKE CARE OF YOUR
LOVED ONES

I - Introduction	4
A - Preamble.....	5
B - General recommendations.....	5
C - Hardware requirements and compatibility.....	6
D - Your Dcare equipment.....	7
II - Your system	10
A - The Dcare ecosystem.....	10
B - Main features.....	11
C - Use of the bracelet.....	16
III - Installation and setup	17
A - Placing the Plugs and learning.....	17
B - Settings and customization.....	20
IV - Miscellaneous	21
A - Indicator LED codes.....	21
B - Replace the battery of the bracelet.....	25
C - FAQ.....	27
D - Warranty.....	28

Thank you for purchasing a Dcare system.

The Dcare solution is a system that ensures the well-being of loved ones and the maintenance of intergenerational ties in families.


Simple, discreet and effective, the Dcare solution ensures good ageing thanks to a connected bracelet.

From a mobile application, the wearer's family and friends are reassured via activity monitoring and are alerted in the event of a problem. Seniors can therefore continue to live at home in total autonomy and peace of mind.

This user manual guides you through the few steps of installation and explains in detail how the solution works, for an optimal experience.

A – Preamble

Because the Dcare system is constantly evolving to keep pace with innovation and your expectations, some of the material in this manual may seem outdated. In this case, up-to-date information is available on our website.

B – General recommendations

Coverage and indoor range

The Dcare solution aims to ensure the well-being of seniors inside their homes. In order to do this, the system needs to cover the entire premises to ensure monitoring and detection at all times.

Thanks to a few smart Plugs placed at various strategic locations inside the house, it ensures the safety of the wearer in every room.

For full coverage, these Plugs should be placed on direct wall sockets that are not obstructed by furniture or substantial obstacles.

The signal range depends on many factors such as the presence of walls and their material. In case the house is not fully covered by the basic set of 3 Dcare Plugs, additional Plugs are available for rent as an option.

C – Hardware requirements and compatibility

Internet connection

In addition to the basic Dcare kit provided at the beginning of the rental contract, the system needs an internet connection in the wearer's home in order to feed information from the Dcare bracelet to the mobile application.

Application

The Dcare system works via an application available for free on Android and iOS systems.

Equipment

Contents of the Dcare Starter Kit:

- A Dcare bracelet
- A Dcare Box with its power supply and RJ45 Ethernet cable
- 3 Dcare Plugs (smart plugs) with integrated 230V classic socket
- A Torx tool
- A spare battery for the bracelet


D – Your Dcare equipment

Dcare Box

Central point of the system, connected to an Internet connection.

It is the relay between the home installation and the mobile application to collect the data captured by the Dcare bracelet.

- Power supply: 5.TV - 2A
- Network connection via Ethernet port

How to install it?

The Dcare Box must communicate with your Dcare Plugs to receive tracking information. Therefore, make sure that it is positioned in an open area, on a piece of furniture for example, close to your Internet router.

Then simply plug it into the mains via the supplied power cable and connect it to the router with the Ethernet cable also supplied.


Understand the Dcare Box

The Dcare Box provides information via light signals, which are detailed in the **Indicator LED codes** section (p.21).

Dcare Plug

Relay between your Box and your bracelet all around the home, they are plugged in wall sockets in strategic spots.

- 2500W (resistive load)
- Input: 230V AC (+/- 15%) 50 Hz
- Directives:
EMC-directive: 2014/30/UE
LV-directive: 2014/35/UE

How to install them?

In order to achieve optimum coverage throughout the home, the 3 (or optional more) Dcare Plugs must be logically placed and configured. For more details, see the **Placing the Plugs and Learning** (p.17).

⚠ Note : The Dcare plugs are also classic 230V electrical sockets. However, for safety reasons, this function is originally disabled. The socket can be switched on by pressing the LED indicator button (which turns red).

Dcare bracelet

It allows the monitoring of the wearer's activity (position and life habits) and ensures safety at all times thanks to its manual alerts and its detection of heavy and soft falls.

- Battery: CR2450 - 3V - 620mAh
- Waterproof
- Low-power


How to configure it?

Before equipping your loved one with the Dcare bracelet, you need to set it up and link it to your system and application. More details in the **Installation guide** (p.19) available on our website.

Understand the Dcare bracelet

The indicator LED on the Dcare bracelet provides you information via light codes. These are detailed in the **Indicator LED codes** section (p.21).

A - The Dcare ecosystem

B – Main features

All the elements of the Dcare solution work together in a complete ecosystem to give you access to a wide range of features. Here is a list of the main features and a description of what they allow.

Manual alerts

In the event of a problem, the wearer can send a manual distress signal via the Dcare bracelet to notify the supervisor(s). To do this, the wearer simply presses and holds the button on the bracelet for two seconds.

Any application connected to this Dcare bracelet will then receive a notification stating that the wearer needs help. Whoever receives this notification can then trigger the necessary actions. For example, the Dcare application offers to call the number registered during the configuration of the installation and corresponding to the wearer (classic GSM call).

Fall detection

The Dcare bracelet is equipped with 5 on-board sensors as well as a smart algorithm, allowing it to detect the wearer's falls. It detects both heavy and soft falls, in order to prevent any type of danger.

When a fall is detected, the wearer can press the main button on the wristband to cancel the alert (e.g. in the case of a fall with recovery and no consequences). If no action is taken within a few seconds of the fall, the information is sent to the Dcare Box and a notification is sent to the supervisor(s)' smartphone(s).

The notification can be treated in the same way as a manual alert but will be qualified differently in the application. The image displayed and the status will indicate a fall of the wearer.

Activity monitoring

Thanks to artificial intelligence (AI), the system analyzes the wearer's daily habits and provides a warning in case of a break in routine via a notification from the application.

A simple example: your parent is still in bed at 11am when he or she is usually up at 8am. You are then directly notified that something is wrong. You can then take any action you wish to remedy the problem. The application offers to call (classic GSM call) the number directly linked to the wearer, registered during the configuration of the installation and the bracelet.

Position request

At any time, a supervisor can request the wearer's last known position through the bracelet data (updated every 15 minutes). To do so, a position icon is available on the application's monitoring screen.


Exit alerts

The smart Plugs installed in the home operate within a limited range. This allows them to operate within the confines of the building and thus cover the whole house efficiently.


When the bracelet is out of range of the installation, the Dcare system considers this as leaving the home. In order to protect the most vulnerable wearers who may be disorientated, the system sends a notification via the application to alert the supervisor(s) that the wearer has left the house. This feature allows frail people to be quickly taken care of when an exit is potentially dangerous.

A notification is also sent to inform of the return of the wearer home after an exit.

Long battery life

For all these features, the Dcare bracelet offers up to 6 months of autonomy. The Dcare Plugs and the Dcare Box are constantly powered from their respective sockets and the Dcare bracelet uses a high-quality battery: no regular recharging is required. A second battery is already included in the Dcare starter kit for initial replacement.

On the monitoring screen of the Dcare app, an indicator shows how much power is left in the battery at any given time. When the battery reaches its end of life, a notification is automatically sent via the app and the indicator will be highlighted. This saves the wearer from worrying about the battery life of the bracelet and helps to anticipate battery changes.

C - Use of the bracelet

In order to configure it and to use its manual alert functionalities, the Dcare bracelet calls on its button according to a few manipulations:

Manual help request : Press and hold the main button until the LED turns red (approximately 2 seconds). It will then continue to flash red until an alert cancellation or after a 3 minutes delay.

Each request for manual assistance and fall detection can be cancelled by the wearer in the event of a false alarm or accidental press within 30 seconds.

Alert cancelation : Press the button twice in succession when the LED is flashing red to cancel an alert (manual help request or fall detection).

Switching to configuration mode: press the button 4 times in a row, holding it down for the 4th time and releasing it when the LED turns blue. It then flashes blue.

This mode is only useful during the configuration of the Dcare bracelet, at the end of which it will exit configuration mode by itself. You should never have to exit the configuration mode manually.

Bracelet reset (in configuration mode): press 6 times in a row to reset the bracelet while in its configuration mode.


This feature is only useful in case of problem during the configuration and should only be used when recommended by the Dcare support.

A – Placing the Plugs and learning

Placing the Plugs

For optimal operation and accurate indoor location, it is imperative to place the Dcare Plugs in the home with the following recommendations in mind:

1. Move the Dcare Plugs as far apart as possible to cover the largest area in the home, in the idea of triangulation.

2. Avoid placing the Dcare Plugs in rooms that are not used regularly.
3. Position the Dcare Plug on an wall, preferably not obstructed by furniture or any other obstacle.

4. Use a direct wall socket rather than a power strip.

There is a risk that a power strip will be moved - and therefore the indoor location will not be accurate - or that it will be switched off by accident. If you want to use a power strip, make sure that it is never moved or turned off.

Also avoid connecting a power strip to the Dcare Plug socket. If the socket has a lot of appliances plugged into it, an overload could damage the plug.

5. Do not move the Dcare Plug after the learning process.

If a Dcare Plug is moved after learning, then it will be necessary to redo the learning process to avoid indoor location problems.

6. If the location is not performing well enough, adding Plugs can enhance the precision.

Learning process

Recommendations

1. Only enter rooms that are regularly used in the list (living room, kitchen, bathroom, etc.). Rooms with passage or where little time is spent (night hall, laundry room, pantry, etc.) should not appear in your list.
2. When moving around the room to determine the boundaries, it is recommended to follow a “snail’s shell” trajectory. Within the time limit, walk along the walls and then gradually move towards the center of the room. In this way you make sure to cover the entire room.

Installation guide

You will also find on our website an **Installation guide** for an installation guided step by step.

B - Settings and customization

In order to optimize battery life and avoid necessary notifications, you can disable certain functions of your system via the Dcare by Domintell app. Some features may not be relevant for your loved one and we allow you to enable only those that are.

As Dcare is an evolving solution, we invite you to visit our website to discover the latest settings and customizations available to optimize your experience with the Dcare system.

A – Indicator LED codes

Dcare bracelet

In order to inform you about the alerts issued or the active operating mode, the Dcare bracelet informs you via its LED indicator:

- **Fast blue flashing:** the bracelet is in configuration mode, ready to be configured.
- **Slow blue flashing:** the bracelet is in configuration mode, it has been provisioned (partially configured).

The configuration mode is only activated by means of a specific operation and is only used when configuring the installation. The handling has been made complex in order to avoid accidental activations.

- **Flashing red:** an alert has been detected and sent (fall or request for help).

It confirms that an alert has been detected and that the supervisor(s) is/are notified. It occurs 30 seconds after the detection (cancellation delay).

- **Very short green flash:** visual aid to identify each button press.

- **Short green flash then blank:** cancellation of alert (request for help or fall).

In case of detection of an alert or sending of a manual alert, the wearer has 30 seconds to decide to cancel the alert by pressing the button on the Dcare bracelet twice.

Dcare Plug

To inform you about its status and its ability to receive information from the Dcare bracelet, the Dcare Plug informs you via its indicator LED:

- **Steady red:** the Dcare Plug is switched on and is correctly configured.

Off: The Dcare Plug is operating and configured but it is not supplying power to its female socket.

- **Slowly flashing red:** The Dcare Plug is not configured.

Dcare Box

To inform you about the alerts issued or the active operating mode, the Dcare Box informs you via its LED indicator:

LED indicator

- **Slowly flashing green:** the Dcare Box is configured and functional. This is the usual status.
- **Steady blue:** the Dcare Box is starting up.

- **Steady yellow:** the Dcare Box cannot access Internet: no network connection, unplugged Ethernet cable or no Cloud service

Also yellow before any QR code has been scanned through the app when the Dcare Box has never been configured.

- **Slowly flashing blue:** the Dcare Box is configured but the learning process has not been completed.

- **Slowly flashing red:** the Dcare Box is not yet configured but is ready for configuration.

- **Rapidly flashing blue:** the configuration of the Dcare Box is incomplete. This may be caused by the app being closed prematurely during a configuration step. In this case, you must restart the configuration in the app.

If errors are repeatedly displayed in your application (especially if a Plug/bracelet was configured), contact us. We will then guide you through a complete system reset.

Advanced

- **Steady red:** there was a problem creating the configuration files, the Dcare Box is not ready for configuration. Unplugging and replugging the Dcare Box will solve the problem.
- **Rapidly flashing red:** a database file is missing. Unplugging and replugging the Dcare Box will solve the problem.

B – Replace the battery of the bracelet

To make changing the battery easier, a Torx tool is included in the Dcare starter kit.

1. Unscrew the 4 screws on the back of the Dcare bracelet case.
2. Remove half of the plastic case and then remove the electronic board by pulling slightly on the metal part containing the battery.
3. Remove the battery from its metal container. Feel free to use the Torx tool provided to push the battery out of its container. Pay attention to the up side of the battery.
4. Replace it with the new battery (CR 2450 - 3V) by inserting it instead of the old one. The side with the '+' sign should be pointing upwards. Do not hesitate to push with a little force and tilt the battery to get it under the contact paddles.

5. Put the circuit board back into the housing, firmly seated at the bottom between the white pins. The hole for the LED and the holes for the two gold contacts (see photo below) are indicators of the way how to set it.

6. Once the electronic card has been replaced, the case can be closed by simply pressing in the direction of the strap closure towards upwards. If it does not close properly directly, the card is incorrectly positioned or in the wrong direction.

7. Tighten the 4 screws on the back of the case.

A short video on our website details this maneuver. Visit dcare.be or scan the QR code on the opposite page.

C – FAQ

We have compiled the most frequently asked questions about our system, its technology and its operation in a **Frequently Asked Questions** section available on our website. Here, you will find answers to any questions or concerns you may have.

dcare.be

As our technology will continue to develop and evolve, we will regularly update our website and this page in particular, in order to answer all your questions.

D – Warranty

The Dcare system and all its products are not medical solutions or health solutions. As a purely informative technology, we decline all responsibility in case of accident and/or physical or moral harm to any person involved in the use of the Dcare system or its services.

All Dcare equipment must be installed and operated in accordance with the manufacturer's recommendations and the standards and regulations for construction and electrical installations in force at the time of installation.

Box and Plugs

Dcare Box and Plug modules are covered by a 2-year warranty.

The Dcare Box must only be powered with the equipment provided. No guarantee will be granted if the Box and Plug modules are transformed, are maintained contrary to the prescribed rules (e.g. open boxes, oxidation, etc.), suffer damage caused by external natural factors (e.g. lightning, frost, water infiltration, humidity, etc.) or damage caused by external accidental factors (e.g. fall, fire, overloading, erroneous command, etc.). The warranty does not apply to any damage resulting from use contrary to Dcare's recommendations or any other cause unrelated to the Dcare software or hardware.

Bracelet

The Dcare bracelet is covered by a 2-year warranty.

No warranty will be granted if the Dcare bracelet is altered, is maintained contrary to the prescribed rules (e.g. oxidation, etc.), suffers damage caused by external natural factors (e.g. lightning, frost, etc.) or damage caused by external accidental factors (e.g. fire, overloading, incorrect control, etc.). The warranty does not apply to any damage resulting from use contrary to Dcare's recommendations or any other cause unrelated to the Dcare software or hardware.

Domintell Cloud

The Dcare system works by combining physical modules with an online data storage, transmission and processing service called «Cloud Domintell». Domintell is in charge of processing and automating the data that is transmitted.

The Cloud Domintell is hosted by a third-party provider (AWS), which ensures the online availability and maintenance of the service. Dcare is not responsible in the event of a failure of the host and its service, which would temporarily interrupt the Cloud service Domintell and therefore the proper functioning of your system. However, we undertake to do our utmost to restore the system as quickly as possible.

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